

HORIBA
Scientific

LAQUA



| | | | |
|-------------|------------------------|------------------|--------------|
| pH | ORP | Ion | Conductivity |
| Resistivity | Total Dissolved Solids | Dissolved Oxygen | Salinity |

ELECTRODES



www.horiba-laqua.com



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

LAQUA Electrode Technology

Born from the fusion of our technical expertise and state-of-the-art manufacturing

As a leading pH electrode manufacturer, HORIBA uses the latest technology for all your measurement needs.

Since the development of Japan's first glass electrode for pH meter, HORIBA has focused on continually improving our electrode technology, especially in materials and manufacturing. HORIBA is committed to continually explore and employ groundbreaking solutions in manufacturing next-generation electrodes so that we always provide you with the newest and best electrodes.

| pH Electrode Selection Guide | | 3-in-1 ELECTRODES | | | | | | | | |
|------------------------------|-----------------------------------|-------------------|----------|----------|----------|----------|----------------|------------|-------------|--------------|
| | | PLASTIC | | | | | STANDARD Tough | LONG Tough | MICRO Tough | SLEEVE Tough |
| | | 9651-10D | 9625-10D | 9630-10D | 9631-10D | 9632-10D | 9615S-10D | 9680S-10D | 9618S-10D | 9681S-10D |
| Specification | Applicable temperature range (°C) | 0-60 | 0-100 | 0-100 | 0-60 | 0-100 | 0-100 | 0-100 | 0-60 | 0-60 |
| | Diameter (mm) | 16 | 16 | 16 | 16 | 16 | 12 | 8 | 3 | 12 |
| | Length (mm) | 150 | 150 | 150 | 155 | 150 | 198 | 283 | 185 | 203 |

| pH - Sample Conditions | | | | | | | | | | |
|------------------------|--|--------------------------------|---|---|---|---|---|---|---|---|
| Aqueous Solution | Conductivity | Normal (over 100 mS/m) | ● | ● | ● | ● | ● | ● | ● | ● |
| | | Low (approx. 10 ~100 mS/m) | | | ● | | | | | ○ |
| | | Very low (approx. 5 ~100 mS/m) | | | ○ | | | | | ○ |
| | | High (approx. 5 S/m) | ○ | ○ | ○ | ○ | ○ | ○ | | ● |
| | Strong alkaline (pH 10-12) | | | | | ● | ○ | ○ | | ○ |
| | Strong acidity (pH 0-2) * Except HF sample | | | | | ● | ● | | | |
| | Quick heat change (within 50°C) | | ● | ● | ● | ● | | | | |
| | High viscosity (approx. 5 Pa·S) | | | | | | | | | ● |
| | Containing non-aqueous solvent | | | | | | ○ | ○ | ○ | ○ |
| | Suspension | | | | | | ○ | ○ | ○ | ● |
| Solid/ Semisolid | Inside | | | | | | | | | |
| | Surface | | | | | | | | | |

| | | | | | | | | | | |
|-------------------|---------------------------|--------------------------|---|---|---|---|---|---|---|---|
| Sample Containers | Microtube/plate (> 50 µL) | | | | | | | | ● | |
| | Ampule | > Ø4 mm | | | | | | | ● | |
| | Micro container (> 2 mL) | | | | | | | ○ | ● | |
| | Tube | ID:13 mm, L:100 ~ 150 mm | | | | | | ● | | |
| | Beaker | 10 mL ~ 1 L | ● | ● | ● | ● | ● | ○ | ○ | ○ |
| | Large container (> 1 L) | | ○ | ○ | ○ | ○ | ○ | ● | | |
| | Petri dish | | | | | | | | | |
| | Droplet | | | | | | | | | |

| | | | | | | | | | | |
|-----------------------------------|--|--|---|---|---|--|---|---|---|---|
| Water | Pure/ion-exchange water (approx. 0.1 mS/m)/ Distilled water (approx. 0.5 mS/m) | | | | | | ○ | | | |
| | Tap/drinking water (approx. 10 mS/m) | | ○ | ○ | ● | | ○ | | | ○ |
| | Surface water | | | | ● | | ○ | | | ○ |
| | Pharmaceutical water/ Environmental water/acid rain | | ○ | ○ | ○ | | ○ | | | ○ |
| Chemical reagent/ solvent | Caustic/strong acid (Except HF sample) | | | | ● | | ● | | | ○ |
| | Hydrofluoric acid | | | | ● | | | | | |
| | Surfactant | | | | | | ○ | | | ● |
| | Water-based paint | | | | | | ○ | | | ● |
| | Dye/coloring agent | | | | | | | | | ● |
| Pharmaceutical/ biological sample | Protein-containing sample | | | | | | ○ | | ○ | ● |
| | Medicinal preparation | | | | | | | | ○ | ○ |
| | Enzyme solution | | | | | | | ○ | ● | |
| | Tris buffer | | | | | | ● | | ○ | ○ |
| | Suspension | | | | | | ○ | | | ● |
| | Agar medium | | | | | | | | | |
| Food | Jam | | | | | | ○ | | | ● |
| | Meat/fish/Fruit/vegetable/ Dough | | | | | | | | | |
| | Honey | | | | | | | | | |
| | Cheese/butter | | | | | | | | | |
| | Yogurt | | ○ | ○ | ○ | | ○ | | | ○ |
| Beverage/ seasoning | Beer | | ○ | ○ | ○ | | ○ | | | ● |
| | Milk/Carbonated drink/juice/ sauce/soy sauce | | | | | | ○ | | | ● |
| | Mayonnaise/ketchup | | | | | | ○ | | | ● |
| Cosmetic/ lotion | Beauty cream/mascara | | | | | | ○ | | | ● |
| | Gel/soap/shampoo/Hairdye lotion | | | | | | ○ | | | ● |
| | Emulsified liquid | | | | | | ○ | | | ○ |



| | | | | COMBINATION ELECTRODES | | | | | | | ISFET ELECTRODES |
|--|----------|-------------|----------|------------------------|----------------|-------------|--------------|----------|------------|----------|---------------------|
| | SLEEVE | NON-AQUEOUS | NEEDLE | PLASTIC | STANDARD ToupH | MICRO ToupH | SLEEVE ToupH | LONG | LONG ToupH | FLAT | GENERAL |
| | 6367-10D | 6377-10D | 6252-10D | 9425-10C | 9415-10C | 9418-10C | 9481-10C | 6069-10C | 9480-10C | 6261-10C | 0040-10D |
| | 0-60 | 0-60 | 0-60 | 0-100 | 0-100 | 0-60 | 0-60 | 0-60 | 0-100 | 0-50 | 0-60 |
| | 12 | 12 | 12 | 16 | 12 | 3 | 12 | 3 | 8 | 12 | 16 |
| | 150 | 150 | 150 | 150 | 198 | 185 | 203 | 291 | 283 | 150 | 190 |

| | | | | | | | | | | |
|---|---|---|---|---|---|---|---|---|---|---|
| ● | ● | ● | ● | ● | ● | ● | ○ | ● | ● | ● |
| | ● | | | | | | ○ | | | |
| | ● | | | | | | ○ | | | |
| | | | ○ | ○ | | ● | | ○ | | |
| ○ | | | | ○ | | ○ | | ○ | | |
| | | | | ● | | | | | | |
| | | | ● | | | | | | | |
| ○ | ● | | | | | ● | | | | |
| ○ | ● | | | ○ | ○ | ○ | | ○ | | ○ |
| | ● | | | ○ | ○ | ● | | ○ | | ○ |
| | | ○ | | | | | | | ● | ● |

| | | | | | | | | | | |
|---|---|---|---|---|---|---|---|---|---|---|
| | | | | | ● | | | | | |
| | | | | | ● | | ○ | | | |
| | | | | | ● | | ○ | ○ | | |
| | | | | | | | ● | ● | | |
| ○ | ○ | ○ | ● | ● | ○ | ○ | ○ | ○ | ○ | ○ |
| | | | ○ | ○ | | | | ● | | |
| | | | | | | | | | ● | ● |
| | | | | | | | | ● | ● | ● |

| | | | | | | | | | | |
|---|---|---|---|---|---|---|--|---|---|-------------|
| | ● | | | ○ | | | | | | |
| | ● | | ○ | ○ | | ○ | | | | |
| | ● | | | ○ | | ○ | | | | |
| | ○ | | ○ | ○ | | ○ | | | | |
| | | | | ● | | ○ | | | | |
| | ○ | | | ○ | | ● | | | | |
| | ○ | | | ○ | | ● | | | | |
| | ○ | | | ○ | | ● | | | | |
| ○ | | | | ○ | ○ | ○ | | | | |
| | ○ | | | | ○ | ○ | | | | |
| | | ○ | | | ● | | | ○ | | |
| | ● | | | ○ | | ○ | | | | |
| | ○ | ○ | | ○ | | ● | | | ● | ● (surface) |
| | | ● | | | | | | ○ | | ● (surface) |
| | ● | | | | | | | | | ○ (surface) |
| | | ○ | | | | | | ○ | | ○ (surface) |
| ○ | | ○ | ○ | ○ | | ○ | | ○ | | ● (surface) |
| ○ | ● | | ○ | ○ | | ● | | | | |
| ○ | ○ | | | ○ | | ● | | | | |
| | ○ | | | ○ | | ● | | | | |
| | ○ | ○ | | ○ | | ● | | | | |
| | ○ | | | ○ | | ● | | | | |
| | ● | | | ○ | | ○ | | | | |



Expertise in Manufacturing

Sophisticated processing technology

HORIBA's in-house expertise in the manufacture of electrodes is the accumulation of more than 60 years of experience. Our sophisticated electrode processing technology provides flexibility in designing various shapes of the electrode bulb and different structural designs of the electrodes.

Thick membrane technology

HORIBA's glass moulding technology allows the manufacture of tougher pH glass bulbs.

Miniaturization



Unique flat electrode design as well as 3mm diameter micro-electrode with integrated temperature sensor (US Patent No. 7314541/ China Patent No. ZL0315796)

Fast response & highly accurate

ToupH glass bulb does not compromise responsiveness and sensitivity (US Patent No. 8262877). Specially designed electrodes are available for hydrofluoric acid & strong alkaline application.

Double-junction electrodes

All HORIBA pH combination electrodes are double-junction electrodes. Flexible to use in a wide-range of applications.

Convenient slider

Refillable electrodes are equipped with a slider to open or close the refilling port easily.

Built-in clip for hooking onto electrode stand arm

Top housing of electrodes is designed with a built-in clip to hook onto HORIBA's electrode stands.

ORP Electrode

| Model | Part No. | Material | Temp. Range (°C) | Application |
|----------|------------|------------|------------------|---|
| 9300-10D | 3014046710 | Pt / Glass | 0 - 60 | Waterproof; Platinum on the flat tip allows measurement of small volume samples |

Ion Selective Electrodes (ISEs)

| Model | Part No. | Combination ISE | Temp. Range (°C) | Measurement Range | Replacement Tip | Part No. |
|-----------|------------|---|------------------|---|-------------------------------|------------|
| 5002S-10C | 3200698386 | Ammonia (NH ₃) | 0 - 50 | 0.1 - 1,000 mg/L NH ₃ | NH ₃ Membrane Caps | 3200705774 |
| 6583S-10C | 3200697410 | Calcium (Ca ²⁺) | 0 - 50 | 0.4 - 40,080 mg/L Ca ²⁺ | 7683S | 3200697414 |
| 6560S-10C | 3200697407 | Chloride (Cl ⁻) | 0 - 50 | 0.35 - 35,000 mg/L Cl ⁻ | 7660S | 3200697411 |
| 6561S-10C | 3200693774 | Fluoride (F ⁻) | 0 - 50 | 0.2 - 19,000 mg/L F ⁻ | 7661S | 3200693606 |
| 6581S-10C | 3200697408 | Nitrate (NO ₃ ⁻) | 0 - 50 | 0.62 - 62,000 mg/L NO ₃ ⁻ | 7681S | 3200697412 |
| 6582S-10C | 3200697409 | Potassium (K ⁺) | 0 - 50 | 0.04 - 39,000 mg/L K ⁺ | 7682S | 3200697413 |

Conductivity Cells

| Type | Model | Part No. | Temp. Range (°C) | Cell Constant | Measurement Range | Application |
|-------------|----------|------------|------------------|----------------------|----------------------|---|
| Submersible | 3551-10D | 3014081712 | 0 - 60 | 0.1 cm ⁻¹ | 0.1 µS/cm - 10 mS/cm | Low conductivity water (e.g., deionized, distilled) |
| | | | | 10 m ⁻¹ | 10 µS/m - 1 S/m | |
| | 9382-10D | 3014046709 | 0 - 80 | 1 cm ⁻¹ | 1 µS/cm - 100 mS/cm | General purpose use; Waterproof |
| | | | | 100 m ⁻¹ | 0.1 mS/m - 10 S/m | |
| | 3552-10D | 3014081545 | 0 - 100 | 1 cm ⁻¹ | 1 µS/cm - 100 mS/cm | General purpose use |
| | | | | 100 m ⁻¹ | 0.1 mS/m - 10 S/m | |
| Flow | 3553-10D | 3014081714 | 0 - 60 | 10 cm ⁻¹ | 10 µS/cm - 1 S/cm | High conductivity water |
| | | | | 1000 m ⁻¹ | 1 mS/m - 100 S/m | |
| | 3561-10D | 3014082350 | 0 - 60 | 0.1 cm ⁻¹ | 0.1 µS/cm - 10 mS/cm | Low conductivity water (e.g., deionized, distilled) |
| | | | | 10 m ⁻¹ | 10 µS/m - 1 S/m | |
| | 3562-10D | 3014082513 | 0 - 60 | 1 cm ⁻¹ | 1 µS/cm - 100 mS/cm | General purpose use |
| | | | | 100 m ⁻¹ | 0.1 mS/m - 10 S/m | |
| | 3573-10C | 3014082590 | 0 - 60 | 10 cm ⁻¹ | 10 µS/cm - 1 S/cm | High conductivity water |
| | | | | 1000 m ⁻¹ | 1 mS/m - 100 S/m | |
| | 3574-10C | 3014082592 | 0 - 60 | 10 cm ⁻¹ | 10 µS/cm - 100 mS/cm | Small volume sample (e.g., column chromatography) |
| | | | | 1000 m ⁻¹ | 1 mS/m - 10 S/m | |

• Material: All have platinum-platinum black / glass-body, except 9382-10D (titanium-platinum black / plastic-body).







Dissolved Oxygen Probes


| Type | Model | Part No. | Temp. Range (°C) | Measurement Range | Replacement Tip | Part No. |
|-------|-----------|------------|------------------|-------------------|-----------------|------------|
| Field | 9551-20D | 3014047090 | 0 - 40 | 0 - 19.99 mg/L DO | 5401 | 3014072770 |
| | 9551-100D | 3014047091 | 0 - 40 | 0 - 19.99 mg/L DO | 5401 | 3014072770 |
| Lab | 9520-10D | 3014046711 | 0 - 45 | 0 - 19.99 mg/L DO | 7541 | 3014074145 |

pH Combination Electrodes (G, R)

HORIBA pH Combination electrodes manufactured with 1 meter cable terminating in BNC connector allow these electrodes to be used with any pH meter¹. Enjoy the full spectrum of features and benefits of these electrodes on your existing pH meter¹. (For applications where temperature measurement and compensation is required, please refer to the 3-in-1 pH electrodes).

¹ pH meters must have BNC connector

| Model | pH Range | Operating Temperature Range (°C) | Liquid Junction | Application |
|---|----------|----------------------------------|-----------------|--|
| ToupH Standard Electrode 9415-10C General laboratory application  <p>Overall length: 198 mm Diameter of probe: 12 mm Connector: BNC</p> <p>3200611623</p> | 0-14 | 0-100 | Ceramic | <p>The electrode offers quick stability and drift reduction.</p> <ul style="list-style-type: none"> Constructed with responsive glass that is 10X stronger than JIS standards The one-touch refilling port slider allows one-hand operation Waterproof, Pb-free glass <p>Perfect for preparing pH buffers and other aqueous test solutions.</p> |
| Standard Plastic Electrode 9425-10C General field application  <p>Overall length: 150 mm Diameter of probe: 16 mm Connector: BNC</p> <p>3200611625</p> | 0-14 | 0-100 | Ceramic | <p>The electrode has plastic body, which is ideal for field measurement.</p> <ul style="list-style-type: none"> Can be submerged up to 1m depth and 30mins (with refilling port closed) Waterproof, Pb-free glass <p>Recommended for field use. For measurement of tap water and drinking water.</p> |
| ToupH Sleeve Electrode 9481-10C High viscosity application  <p>Overall length: 203 mm Diameter of probe: 12 mm Connector: BNC</p> <p>3200611631</p> | 0-14 | 0-60 | Movable sleeve | <p>The electrode gives stable readings in highly viscous samples.</p> <ul style="list-style-type: none"> The liquid junction is designed with a movable sleeve that can be cleaned easily and prevents clogging Waterproof, Pb-free glass <p>For measurement of highly viscous samples and samples containing non-aqueous solvents (e.g., cosmetics, paints).</p> |
| ToupH Micro Electrode 9418-10C Precious trace amount sample  <p>Overall length: 185 mm Diameter of probe: 3 mm Connector: BNC</p> <p>3200611627</p> | 0-14 | 0-60 | Ceramic | <p>The electrode can measure samples as small as 50µL.</p> <ul style="list-style-type: none"> Compatible with extremely small containers (e.g., micro tubes) Temperature sensor is placed next to the bulb for quick response Waterproof <p>Suitable for low-volume samples and wide range of aqueous solutions.</p> |
| ToupH Long Electrode 9480-10C For large containers and long test tubes  <p>Overall length: 283 mm Diameter of probe: 8 mm Connector: BNC</p> <p>3200611628</p> | 0-14 | 0-100 | Ceramic | <p>The long, thin body of the electrode is perfect for large containers and test tubes.</p> <ul style="list-style-type: none"> 283mm length, 8mm diameter Constructed with responsive glass that is 10X stronger than JIS standards Waterproof, Pb-free glass <p>For measuring samples (e.g., microbial culture fluids) in test tubes and tall beakers.</p> |
| Long Electrode 6069-10C For very slender test tubes  <p>Overall length: 291 mm Diameter of probe: 3.15 mm Connector: BNC</p> <p>3014081107</p> | 0-14 | 0-60 | Ceramic | <p>The long, thin body of the electrode is perfect for very slender test tubes.</p> <ul style="list-style-type: none"> 291mm length, 3mm diameter Waterproof <p>For measuring samples in slender tubes (e.g., NMR test tube).</p> |






| Model | pH Range | Operating Temperature Range (°C) | Liquid Junction | Application |
|---|----------|----------------------------------|-----------------|--|
| Flat Electrode (G, R) 6261-10C  <p>Overall length: 150 mm Diameter of probe: 12 mm Connector: BNC</p> | 0-12 | 0-50 | Sleeve | <p>The sensor is located on the flat surface of the tip.</p> <ul style="list-style-type: none"> Measurement can be made from minute amount of moisture on solid sample surface Pure water can be applied for samples with no moisture Waterproof <p>Perfect for measuring samples in shallow containers (e.g., petri dishes) and gelatinous materials (e.g., nutrient agar). For surface measurement of meat, paper, skin, and cloth.</p> |

3014081807



3-in-1 pH Glass Body Electrodes² (GRT)

HORIBA pH Combination electrodes with an integrated thermistor offer higher accuracy as these electrodes measure temperature concurrently with pH. The pH meter is able to continuously monitor and compensate for temperature effects automatically.

² Only compatible with HORIBA pH meters






| Model | pH Range | Operating Temperature Range (°C) | Liquid Junction | Application |
|---|----------|----------------------------------|-----------------|---|
| ToupH Standard Electrode 9615S-10D General laboratory application  <p>Overall length: 198 mm Diameter of probe: 12 mm Connectors: BNC & phono jack</p> | 0-14 | 0-100 | Ceramic | <p>The electrode offers quick stability and drift reduction.</p> <ul style="list-style-type: none"> Constructed with responsive glass that is 10x stronger than JIS standards The one-touch refilling port slider allows one-hand operation Waterproof, Pb-free glass <p>Perfect for preparing pH buffers and other aqueous test solutions.</p> |
| ToupH Sleeve Electrode 9681S-10D High viscosity application  <p>Overall length: 203 mm Diameter of probe: 12 mm Connectors: BNC & phono jack</p> | 0-14 | 0-60 | Movable sleeve | <p>The electrode gives stable readings in highly viscous samples.</p> <ul style="list-style-type: none"> The liquid junction is designed with a movable sleeve that can be cleaned easily and prevents clogging Waterproof, Pb-free glass <p>For measurement of highly viscous samples and samples containing non-aqueous solvents (e.g. cosmetics, paints).</p> |
| ToupH Micro Electrode 9618S-10D Precious trace amount sample  <p>Overall length: 185 mm Diameter of probe: 3 mm Connectors: BNC & phono jack</p> | 0-14 | 0-60 | Ceramic | <p>The electrode can measure samples as small as 50µL.</p> <ul style="list-style-type: none"> Compatible with extremely small containers (e.g. micro tubes) Temperature sensor is placed next to the bulb for quick response Waterproof <p>Suitable for low-volume samples and a wide range of aqueous solutions.</p> |
| ToupH Long Electrode 9680S-10D For large containers and long test tubes  <p>Overall length: 283 mm Diameter of probe: 8 mm Connectors: BNC & phono jack</p> | 0-14 | 0-100 | Ceramic | <p>The long, thin body of the electrode is perfect for large containers and test tubes.</p> <ul style="list-style-type: none"> 283mm length, 8mm diameter Constructed with responsive glass that is 10x stronger than JIS standards Waterproof, Pb-free glass <p>For measuring samples (e.g. microbial culture fluids) in test tubes and tall beakers.</p> |
| Needle Electrode 6252-10D For food application  <p>Overall length: 150 mm Diameter of probe: 12 mm Connectors: BNC & phono jack</p> | 0-12 | 0-60 | Ceramic | <p>Needle electrode allows measurement of food samples and aqueous solutions.</p> |

3014080850

| Model | pH Range | Operating Temperature Range (°C) | Liquid Junction | Application |
|---|----------|----------------------------------|-----------------|---|
| Low-Conductivity Electrode 6377-10D For pure water & non-aqueous solvents  Overall length: 150 mm Diameter of probe: 12 mm Connectors: BNC & phono jack 3014093085 | 0-14 | 0-60 | Movable sleeve | Uses a glass membrane that is highly sensitive to low-conductivity water and non-aqueous solvents. |
| Standard Sleeve Electrode 6367-10D  Overall length: 150 mm Diameter of probe: 12 mm Connectors: BNC & phono jack 3014079136 | 0-14 | 0-60 | Sleeve | Uses a sleeve at the liquid junction for improved stability and repeatability. For measuring pH at high accuracy. |

3-in-1 pH Plastic Body Electrodes² (GRT)

² Only compatible with HORIBA pH meters

| Model | pH Range | Operating Temperature Range (°C) | Liquid Junction | Applications |
|---|----------|----------------------------------|------------------------------|---|
| Gel-filled pH Electrode 9651-10D  Overall length: 150 mm Diameter of probe: 16 mm Connectors: BNC & phono jack 3200642020 | 0-14 | 0-60 | Porous sintered polyethylene | The plastic body of the electrode is filled with gel electrolyte. Less maintenance is needed as refilling is not required. <ul style="list-style-type: none"> Can be submerged up to 1m depth of water for 30mins. Waterproof, Pb-free glass Recommended for field use. |
| Standard Plastic Electrode 9625-10D; 9625-20D; 9625-30D For Field  Overall length: 150 mm Diameter of probe: 16 mm Connectors: BNC & phono jack 3200360505; 3200393025; 3200393026 | 0-14 | 0-100 | Ceramic | The electrode has a plastic body which is ideal for field measurement. <ul style="list-style-type: none"> Can be submerged up to 1m depth of water for 30mins. (with refilling port closed) Waterproof, Pb-free glass Recommended for field use. For measurement of tap water and drinking water. |
| Hydrofluoric Acid Resistant Electrode 9631-10D  Overall length: 155 mm Diameter of probe: 16 mm Connectors: BNC & phono jack 3200524119 | 2-12 | 0-60 | Ceramic | The electrode can measure 1% hydrofluoric acid solution (at 25°C, immersed at 1min.) for about 1000 times. <ul style="list-style-type: none"> Rolled glass design for long-term reliable measurement and easy maintenance Compliant with Japan's Measurement Act Certification Waterproof, Pb-free glass Suitable for drain water measurement after etching process. |
| Strong Alkali Resistant Electrode 9632-10D  Overall length: 150 mm Diameter of probe: 16 mm Connectors: BNC & phono jack 3200524120 | 0-14 | 0-100 | Ceramic | The alkali-resistant glass membrane has higher resistance and longer stability (about 5X in 0.1mol/L sodium at 60°C, pH 13) than conventional electrodes. <ul style="list-style-type: none"> Waterproof, Pb-free glass Suitable for strong alkali samples such as plating solutions. |
| Standard Plastic Electrode 9630-10D For tap water  Overall length: 150 mm Diameter of probe: 16 mm Connectors: BNC & phono jack 3200528726 | 0-14 | 0-100 | Ceramic | The electrode can measure samples with low conductivity or buffering capacity. <ul style="list-style-type: none"> Made of high purity multicomponent lithium series glass Waterproof, Pb-free glass Suitable for tap water measurement and quality control in water purification plant. Recommended to use with cleaning solution 230. |

Next-Generation Electrode Technology


Semiconductor technology without glass

HORIBA started researching **ISFET (Ion Sensitive Field Effect Transistor)** using semiconductor technology many years ago and continued to improve its quality. This has provided a new solution for environments where glass material cannot be used.



ISFET

Applicable electrode: 0040-10D



What is an ISFET (semiconductor sensor)?

ISFET is the abbreviation of **Ion Sensitive Field Effect Transistor**.
The response part uses a semiconductor based sensor.

Special features of the ISFET

1. Will not crack or break like conventional glass electrodes
2. The sensor is flat and very small enabling the measurement of extremely small samples
3. Easy handling and maintenance - simply clean with a toothbrush
4. Can be stored dry

The flat electrode has a distance of less than 100 μm between the housing and sensor


The unique structure allows measurements to be taken from the smallest amount of moisture on solid objects and prevents bubbles being trapped on the sensor when measuring samples in a beaker.

Reduction of static electricity effect


The combination of HORIBA's unique semiconductor device structure together with the improved electrostatic protection circuit results in a significant reduction of the static electricity effect that had previously been the weak point of a semiconductor sensor.

ISFET ELECTRODES **ISFET**

ISFET is the abbreviation of Ion Sensitive Field Effect Transistor. Since ISFET is robust and will not crack like the conventional glass electrodes, it can be easily handled and maintained. The response part is equipped with a flat and miniature semiconductor-based sensor, which makes the measurement even on extremely small samples possible. Combination of HORIBA's unique semiconductor device structure and improvement of the electrostatic protection circuit enables to reduce greatly the static electricity effect that had been the weak point of the semiconductor sensor. Now the measurement has become more comfortable and reliable.

| Model | pH Range | Operating Temperature Range (°C) | Liquid Junction | Applications |
|---|----------|----------------------------------|------------------------------|--|
| General ISFET pH electrode 0040-10D  <small>Overall length: 190.6 mm Diameter of probe: 16 mm Connectors: BNC & phono jack</small> <small>3200367925</small> | 0-14 | 0-60 | Porous sintered polyethylene | <p>The sensor is located on the flat surface of the tip (<100μm from the housing).</p> <ul style="list-style-type: none"> • Measurement can be made from minute amount of moisture on solid sample surface • Pure water can be applied for samples with no moisture • Use of semiconductor sensor prevents damage such as crack or breakage • Waterproof • Replacement sensor (0141) 3200367926 <p>Perfect for measuring samples in shallow containers (e.g., petri dishes) and gelatinous materials (e.g., nutrient agar). For surface measurement of meat, paper, skin, and cloth.</p> |

Metallic Electrode (For ORP Measurement)





| Model | Operating Temperature Range (°C) | Electrode Material | Internal Solution | Applications |
|--|----------------------------------|--------------------|-------------------|---|
| ORP Electrode 9300-10D Waterproof platinum 3-in-1 type  <small>Overall length: 150 mm Diameter of probe: 12 mm Connectors: BNC & phono jack</small> <small>3014046710</small> | 0-60 | Pt / Glass | #300 (KCl) | Waterproof; Platinum on the flat tip allows measurement of small volume samples |

Conductivity Electrode Cells





HORIBA Conductivity cells are available as Submersible type and Flow type, as well as in a variety of cell constants ranging from 0.1 to 10.0.

The HORIBA Conductivity cells are integrated with temperature sensor (except for 3573 & 3574) and the wetted material is **Platinum / Titanium, coated with Platinum black**. Rugged Titanium allows cell to be used in a wide range of applications, including highly corrosive samples such as concentrated acids and sea water. Maintenance is simple – soak in deionized/demineralized water or with the conditioning solution.

Conductivity Cells (Submersible Type)

| Model | Cell Constant | Measurement Range | Temp. Range (°C) | Cell Material | Thermistor | Minimum Sample Volume (ml) | Application |
|---|----------------------|----------------------|------------------|-----------------------|------------|----------------------------|---|
| 3551-10D  Overall length: 175 mm Diameter of probe: 23 mm Connectors: BNC & phono jack 3014081712 | 0.1 cm ⁻¹ | 0.1 µS/cm - 10 mS/cm | 0 - 60 | Pt-Pt black / Glass | Built-in | 50 | Low conductivity water (e.g., deionized, distilled) |
| | 10 m ⁻¹ | 10 µS/m - 1 S/m | | | | | |
| 3552-10D  Overall length: 150 mm Diameter of probe: 12 mm Connectors: BNC & phono jack 3014081545 | 1 cm ⁻¹ | 1 µS/cm - 100 mS/cm | 0 - 100 | Pt-Pt black / Glass | Built-in | 15 | General purpose use |
| | 100 m ⁻¹ | 0.1 mS/m - 10 S/m | | | | | |
| 3553-10D  Overall length: 175 mm Width of probe: 28 mm Connectors: BNC & phono jack 3014081714 | 10 cm ⁻¹ | 10 µS/cm - 1 S/cm | 0 - 60 | Pt-Pt black / Glass | Built-in | 50 | High conductivity water |
| | 1000 m ⁻¹ | 1 mS/m - 100 S/m | | | | | |
| 9382-10D  Overall length: 150 mm Diameter of probe: 16 mm Connectors: BNC & phono jack 3014046709 | 1 cm ⁻¹ | 1 µS/cm - 100 mS/cm | 0 - 80 | Ti-Pt black / Plastic | Built-in | 20-30 | General purpose use; Waterproof |
| | 100 m ⁻¹ | 0.1 mS/m - 10 S/m | | | | | |







Conductivity Cells (Flow Type)







| Model | Cell Constant | Measurement Range | Temp. Range (°C) | Cell Material | Thermistor | Minimum Sample Volume (ml) | Application |
|--|----------------------|----------------------|------------------|---------------------|------------|----------------------------|---|
| 3561-10D  Overall length: 143 mm Diameter of probe: 18 mm Connectors: BNC & phono jack 3014082350 | 0.1 cm ⁻¹ | 0.1 µS/cm - 10 mS/cm | 0 - 60 | Pt-Pt black / Glass | Built-in | 10 | Low conductivity water (e.g., deionized, distilled) |
| | 10 m ⁻¹ | 10 µS/m - 1 S/m | | | | | |
| 3562-10D  Overall length: 205 mm Diameter of probe: 18 mm Connectors: BNC & phono jack 3014082350 | 1 cm ⁻¹ | 1 µS/cm - 100 mS/cm | 0 - 60 | Pt-Pt black / Glass | Built-in | 16 | General purpose use |
| | 100 m ⁻¹ | 0.1 mS/m - 10 S/m | | | | | |
| 3573-10C  Overall length: 222 mm Diameter of probe: 18 mm Connector: BNC 3014082590 | 10 cm ⁻¹ | 10 µS/cm - 1 S/cm | 0 - 60 | Pt-Pt black / Glass | — | 4 | High conductivity water |
| | 1000 m ⁻¹ | 1 mS/m - 100 S/m | | | | | |
| 3574-10C  Overall length: 136 mm Diameter of probe: 66 mm Connector: BNC 3014082592 | 10 cm ⁻¹ | 10 µS/cm - 100 mS/cm | 0 - 60 | Pt-Pt black / Glass | — | 0.25 | Small volume sample (e.g., column chromatography) |
| | 1000 m ⁻¹ | 1 mS/m - 10 S/m | | | | | |

Combination ISE

Ion-selective electrodes are responsive to concentration of particular ions in the test liquid and are variable-potential electrodes. They are used in conjunction with reference electrodes to measure the concentration of particular ions. HORIBA's years of experience and know-how in this field are behind the wide range of ion electrodes we offer.

When measurements are made using an ion meter, calibrating it with various standard solutions will give direct readings of the ion concentration. Note that since volume-detection level changes with temperature, measurements must be taken at a fixed temperature.

| Model | Accessories Included | Temp. Range (°C) | Measurement Range | pH Range |
|--|---|------------------|--|---|
|  <p>Ammonia ion (NH₃) electrode 5002S-10C 3200698386 Overall length: 161 mm Diameter of probe: 15 mm Connector: BNC</p> | <ul style="list-style-type: none"> • membrane cap, 3pcs • 1000mg/L ammonium ion standard solution, 50ml • 100mg/L ammonium ion standard solution, 50ml • ammonia electrode filling solution, 50ml • ammonia ionic strength adjustor, 50ml • syringe • dropper • protective pipe • manual | 0 - 50 | 0.1 - 1,000 mg/L NH ₃ | Adjust more than pH 12 |
|  <p>Calcium ion (Ca²⁺) electrode 6583S-10C 3200697410 Overall length: 150 mm Diameter of probe: 16 mm Connector: BNC</p> | <ul style="list-style-type: none"> • calcium electrode tip, 2pcs • 1000mg/L calcium ion standard solution, 50ml • 100mg/L calcium ion standard solution, 50ml • calcium electrode filling solution, 50ml • calcium ionic strength adjustor, 50ml • syringe • dropper • protective pipe • manual | 0 - 50 | 0.4 - 40,080 mg/L Ca ²⁺ (10 ⁻⁵ to 1 mol/L Ca ²⁺) | 4.0 mg/L (10 ⁻⁴ mol/L) Ca ²⁺ , pH 5 to 11 |
|  <p>Chloride ion (Cl⁻) electrode 6560S-10C 3200697407 Overall length: 150 mm Diameter of probe: 16 mm Connector: BNC</p> | <ul style="list-style-type: none"> • chloride electrode tip • 1000mg/L chloride ion standard solution, 50ml • 100mg/L chloride ion standard solution, 50ml • chloride electrode filling solution, 50ml • chloride ionic strength adjustor, 50ml • syringe • dropper • protective pipe • water-resistant abrasive sheet • manual | 0 - 50 | 0.35 - 35,000 mg/L Cl ⁻ (10 ⁻⁵ to 1 mol/L Cl ⁻) | 350 mg/L (10 ⁻² mol/L) Cl ⁻ , pH 3 to 11 |
|  <p>Fluoride ion (F⁻) electrode 6561S-10C 3200693774 Overall length: 150 mm Diameter of probe: 16 mm Connector: BNC</p> | <ul style="list-style-type: none"> • fluoride electrode tip • 1000mg/L fluoride ion standard solution, 50ml • 100mg/L fluoride ion standard solution, 50ml • fluoride electrode filling solution, 50ml • fluoride ionic strength adjustor, 50ml • syringe • dropper • protective pipe • manual | 0 - 50 | 0.2 - 19,000 mg/L F ⁻ (10 ⁻⁶ to 1 mol/L F ⁻) | 20 mg/L (10 ⁻³ mol/L) F ⁻ , pH 4 to 10 |
|  <p>Nitrate ion (NO₃⁻) electrode 6581S-10C 3200697408 Overall length: 150 mm Diameter of probe: 16 mm Connector: BNC</p> | <ul style="list-style-type: none"> • nitrate electrode tip, 2pcs • 1000mg/L nitrate ion standard solution, 50ml • 100mg/L nitrate ion standard solution, 50ml • nitrate electrode filling solution, 50ml • nitrate ionic strength adjustor, 50ml • syringe • dropper • protective pipe • manual | 0 - 50 | 0.62 - 62,000 mg/L NO ₃ ⁻ (10 ⁻⁵ to 1 mol/L NO ₃ ⁻) | 62 mg/L (10 ⁻³ mol/L) NO ₃ ⁻ , pH 3 to 7 |
|  <p>Potassium ion (K⁺) electrode 6582S-10C 3200697409 Overall length: 150 mm Diameter of probe: 16 mm Connector: BNC</p> | <ul style="list-style-type: none"> • potassium electrode tip, 2pcs • 1000mg/L potassium ion standard solution, 50ml • 100mg/L potassium ion standard solution, 50ml • potassium electrode filling solution, 50ml • potassium ionic strength adjustor, 50ml • syringe • dropper • protective pipe • manual | 0 - 50 | 0.04 - 39,000 mg/L K ⁺ (10 ⁻⁶ to 1 mol/L K ⁺) | 3.9 mg/L (10 ⁻⁴ mol/L) K ⁺ , pH 5 to 11 |

| Selection Coefficient | Replacement Tip | Electrode Filling Solution | 100mg/L Standard Solution | 1000mg/L Standard Solution | Ionic Strength Adjustor | Applications |
|--|---|----------------------------|---------------------------|----------------------------|---------------------------|---|
| — |  <p>NH₃ electrode membrane caps 3200705774</p> | 500-NH3-IFS 3200697173 | 500-NH4-SL 3200697172 | 500-NH4-SH 3200697171 | 500-NH3-ISA 3200697174 | Agriculture, Soil, Power Station Water, Fish Tanks, Sea Water, Waste Water, Plating Baths, Air / Stack Gases and Biological Cultures or Samples |
| $\text{Fe}^{3+} = 0.1$, Fe^{2+} , $\text{Zn}^{2+} = 1$, $\text{Sr}^{2+} = 50$ Ni^{2+} , $\text{Cu}^{2+} = 70$, $\text{Co}^{2+} = 350$ $\text{Mn}^{2+} = 500$, $\text{Mg}^{2+} = 1,000$ Na^+ , K^+ , Ba^{2+} , $\text{NH}_4^+ = \text{over } 1,000$ |  <p>7683S 3200697414</p> <p>Calcium</p> | 500-CA-IFS 3200697177 | 500-CA-SL 3200697176 | 500-CA-SH 3200697175 | 500-CA-ISA 3200697178 | Agriculture / Plant Tissue, Soil, Water Softening Systems, Boiler Feed Water, Drinking / Mineral Water, Biological Cultures, Dental / Clinical Analysis and Dairy / Food / Beverages Applications |
| $\text{S}_2\text{O}_3^{2-}$, S^{2-} , I^- , Ag^+ , $\text{Hg}^{2+} = \text{Not acceptable}$ $\text{SCN}^- = 0.3$, $\text{MnO}_4^- = 0.1$ $\text{Br}^- = 0.03$ NO_3^- , F^- , HCO_3^- , SO_4^{2-} , $\text{PO}_4^{2-} = 1,000$ |  <p>7660S 3200697411</p> <p>Chloride</p> | 500-CL-IFS 3200697169 | 500-CL-SL 3200697168 | 500-CL-SH 3200697167 | 500-CL-ISA 3200697170 | Agriculture, River / Tap Water, Plant Tissue, Soils, Boiler Feed Water, Clinical Analysis, Sweat, Urine, Cement, Plating Baths and Dairy / Food / Beverages Samples |
| Possible interference when multiply-charged ion (ex. Al^{3+} , Fe^{3+}) coexisted and foamed the complex. |  <p>7661S 3200693606</p> <p>Fluoride</p> | 500-F-IFS 3200697165 | 500-F-SL 3200697164 | 500-F-SH 3200697163 | 500-F-TISAB 3200697166 | Dental / Toothpaste / Mouth Wash, Drinking / Seawater, Wastewater, Air / Stack Gases, Acids, Soils, Food, Biological Fluids, Plant Tissue, Coal, Carbonated Beverages and Bone |
| ClO_4^- , $\text{I}^- = \text{Not acceptable}$, $\text{Br}^- = 2$ $\text{NO}_2^- = 3$, $\text{Cl}^- = 300$ HCO_3^- , H_2PO_4^- , $\text{SO}_4^{2-} = \text{over } 1000$ |  <p>7681S 3200697412</p> <p>Nitrate</p> | 500-NO3-IFS 3200697181 | 500-NO3-SL 3200697180 | 500-NO3-SH 3200697179 | 500-NO3-ISA 3200697182 | Agriculture / Plant Tissue / Fertilizers, Surface / Seawater / Drinking Water, Sewage Effluent, Soils, Meats, Vegetables, Foods / Beverages |
| $\text{Rb}^+ = 0.4$, $\text{Cs}^+ = 3$, $\text{NH}_4^+ = 70$ Li^+ , Na^+ , Mg^{2+} , Ca^{2+} , Sr^{2+} , $\text{Ba}^{2+} = \text{over } 1,000$ |  <p>7682S 3200697413</p> <p>Potassium</p> | 500-K-IFS 3200697185 | 500-K-SL 3200697184 | 500-K-SH 3200697183 | 500-K-ISA 3200697186 | Agriculture / Plant Tissue, Soils, Wastewater, River / Tap Water, Clinical Analysis, Saliva, Serum, Fertilizers, Soils and Wines, Dairy / Foods / Beverages |

Note: Detailed information on standard solutions, ISAs, and filling solutions can be found on page 14




| Ion Electrodes | | *All ion electrodes (except combination electrodes) require a sensor holder for attaching to the electrode stand. *Please be aware of the hindering ion and pH range interference of ion electrodes. *D-73 connects combination type ion electrodes only. | | | |
|---------------------------|-----------|--|--------------------------------|---|------------|
| Electrode Name | Model | Measuring Range | Applicable reference electrode | Interfering ion influence ¹ | Part No. |
| Sodium ion electrode | 1512A-10C | 2.3~230,000 mg/L Na ⁺ | 2565A | K ⁺ , Li ⁺ =10 NH ₄ ⁺ =20 Ca ²⁺ =500 | 3014068526 |
| Cyanide ion electrode | 8001-10C | 0.03~2,600 mg/L CN ⁻ | 2060A • 2565A | S ²⁻ , MnO ₄ ⁻ =N/A I ⁻ =0.1 S ₂ O ₃ ²⁻ =1 | 3014094393 |
| Chloride ion electrode | 8002-10c | 0.4~35,000 mg/L Cl ⁻ | 2565A | S ₂ O ₃ ²⁻ , S ²⁻ , I ⁻ , Ag ⁺ , Hg ²⁺ =N/A SCN ⁻ =0.3 MnO ₄ ⁻ =0.1 Br ⁻ =0.03 NO ₂ ⁻ , F ⁻ , HCO ₃ ⁻ , SO ₄ ²⁻ , PO ₄ ³⁻ =1,000 | 3014094394 |
| Sulfide ion electrode | 8003-10C | 0.3~32,000 mg/L S ²⁻ | 2060A • 2565A | CN ⁻ =N/A S ₂ O ₃ ²⁻ =10 I ⁻ , F ⁻ , Cl ⁻ , PO ₄ ³⁻ , SO ₄ ²⁻ =1,000 | 3014094395 |
| Iodide ion electrode | 8004-10C | 0.01~13,000 mg/L I ⁻ | 2060A • 2565A | MnO ₄ ⁻ , S ²⁻ , CN ⁻ =N/A S ₂ O ₃ ²⁻ =10 NO ₂ ⁻ =100 Br ⁻ =1,000 | 3014094396 |
| Bromide ion electrode | 8005-10C | 0.8~80,000 mg/L Br ⁻ | 2565A | S ₂ O ₃ ²⁻ , I ⁻ , S ²⁻ , CN ⁻ =N/A MnO ₄ ⁻ =1 Cl ⁻ , PO ₄ ³⁻ =100 F ⁻ , NO ₃ ⁻ , SO ₄ ²⁻ =1,000 | 3014094397 |
| Copper ion electrode | 8006-10C | 0.06~6,400 mg/L Cu ²⁺ | 2565A | Fe ²⁺ =0.1 Ni ²⁺ , Na ⁺ =1,000 | 3014094398 |
| Cadmium ion electrode | 8007-10C | 0.1~11,000 mg/L Cd ²⁺ | 2060A • 2565A | Cu ²⁺ , Hg ²⁺ , Ag ⁺ =N/A Pb ²⁺ =0.1 Fe ³⁺ =1 Cr ³⁺ Fe ²⁺ =100 Ni ²⁺ =1,000 | 3014094399 |
| Lead ion electrode | 8008-10C | 2~20,000 mg/L Pb ²⁺ | 2565A | Cu ²⁺ , Hg ²⁺ , S ²⁻ , Ag ⁺ =N/A Fe ³⁺ =0.01 Cr ³⁺ =1 Cd ²⁺ =10 Ni ²⁺ , Mg ²⁺ , Zn ²⁺ =100 NH ₄ ⁺ , K ⁺ =1,000 | 3014094400 |
| Thiocyanate ion electrode | 8009-10C | 0.6~5,800 mg/L SCN ⁻ | 2565A | CN ⁻ , I ⁻ , S ²⁻ , S ₂ O ₃ ²⁻ =N/A Br ⁻ =1 Cl ⁻ =100 | 3014094401 |
| Fluoride ion electrode | 8010-10C | 0.02~19,000 mg/L F ⁻ | 2060A • 2565A | Possible interference when multiply-charged ion (ex. Al ³⁺ , Fe ³⁺) coexisted and foamed the complex. | 3014094439 |
| Silver ion electrode | 8011-10C | 0.01~110,000 mg/L Ag ⁺ | 2565A | Hg ²⁺ =N/A Cu ²⁺ , Cd ²⁺ , Pb ²⁺ , Zn ²⁺ , Mg ²⁺ , Ca ²⁺ , Na ⁺ , K ⁺ =Over 1,000 | 3014094402 |
| Nitrate ion electrode | 8201-10C | 0.62~62,000 mg/L NO ₃ ⁻ | 2565A | ClO ₄ ⁻ =0.03 I ⁻ =0.1 Br ⁻ =2 NO ₂ ⁻ =3 Cl ⁻ =40 F ⁻ =200 CH ₃ COO ⁻ =300 SO ₄ ²⁻ =Over 1,000 | 3014094403 |
| Potassium ion electrode | 8202-10C | 0.04~39,000 mg/L K ⁺ | 2565A | Rb ⁺ =0.4 Cs ⁺ =3 NH ₄ ⁺ =70 Li ⁺ , Na ⁺ , Mg ²⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ =Over 1,000 | 3014094404 |
| Calcium ion electrode | 8203-10C | 0.04~40,080 mg/L Ca ²⁺ | 2060A • 2565A | Fe ³⁺ =0.1 Fe ²⁺ , Zn ²⁺ =1 Sr ²⁺ =50 Ni ²⁺ Cu ²⁺ =70 Co ²⁺ =350 Mn ²⁺ =500 Mg ²⁺ =1,000 Na ⁺ , K ⁺ , Ba ²⁺ , NH ₄ ⁺ =Over 1,000 | 30140688 |

DISSOLVED OXYGEN (DO) ELECTRODE & TIPS



HORIBA Dissolved Oxygen (DO) electrodes are galvanic probes with integrated temperature sensors. With galvanic DO probes, calibration can be performed immediately and in air. The HORIBA DO probes use unique and innovative tips which are replaceable. No need to replace membranes or refill electrolytes.

Two models are available: a Laboratory model (9520) that can be used for BOD measurements, and a Field immersible model (9551) housed in a rugged casing available in 2m and 10m cable configurations. The Laboratory 9520 DO probe is fitted with a rotor as well as an adaptor to facilitate BOD measurements.

Dissolved Oxygen Electrodes

| Model | Measurement Range | Response Time | Temperature Range (°C) | Features |
|--|-------------------|---|------------------------|---|
| 9520-10D For laboratories  Overall length: 184 mm Diameter of probe: 15 mm Connectors: BNC & phono jack 3014046711 | 0-19.99mg/L DO | 20 seconds (90% response time at constant temperature) | 0-45 | Waterproof; It operates with the built-in temperature sensor and replaceable DO tip 7541. |
| 9551-20D For field immersible type (2 m cable)  Overall length: 165 mm Diameter of probe: 32 mm Connectors: BNC & phono jack 3014047090 | 0-19.99mg/L DO | 30 seconds (90% response time at constant temperature) | 0-40 | Waterproof; It operates with the built-in temperature sensor and replaceable DO tip 5401. |
| 9551-100D For field immersible type (10 m cable)  Overall length: 165 mm Diameter of probe: 32 mm Connectors: BNC & phono jack 3014047091 | 0-19.99mg/L DO | 30 seconds (90% response time at constant temperature) | 0-40 | Waterproof; It operates with the built-in temperature sensor and replaceable DO tip 5401. |

Dissolved Oxygen Electrode Tips

| Model | Description |
|--|---|
| 5401  3014072770 | Replacement DO tip for 9551-20D and 9551-100D |
| 7541  Overall length: 26.5 mm Diameter: 15 mm 3014074145 | Replacement DO tip for 9520-10D |

pH Electrode Care and Maintenance Procedures

Your pH electrode will eventually reach the end of its useful life as its performance naturally degrades over time. To maximize the performance of your pH electrode and extend its life span, proper care and regular maintenance are equally required.



Part no. 3014028653
Cleaning Solution 220 - contains 10% thiourea and 1% hydrochloric acid (HCl) for removing inorganic residues on glass membrane and junction



Part no. 3200366771
Cleaning Solution 250 - contains less than 0.5% enzyme protease, less than 0.1% sodium azide, and other ingredients (See SDS) for removing protein residues on glass membrane and junction



Mild detergent



Soft lint-free tissue



Clean water (e.g., tap, distilled or deionized water) in a squirt bottle



Part no. 3200043640
Model 300
3.33 KCl
Internal solution for electrode maintenance.

Refer to the safety data sheet (SDS) of the chemical solution to be used in cleaning and wear the appropriate personal protective equipment for safe handling. Download the SDSs of HORIBA solutions at www.horiba-laqua.com.

Refilling

The pH electrode may be filled with either an ionic liquid solution (refillable or liquid-filled pH electrode) or ionic gel solution (sealed or gel-filled pH electrode). Gel-filled pH electrodes do not require routine refilling and typically require less maintenance than liquid-filled electrodes. Liquid-filled pH electrodes are constructed with refilling port, which is securely covered with a slider. The refilling port allows you to fill or empty the reference chamber.

- To top up or re-fill the reference chamber of liquid-filled pH electrode, push the slider upward to uncover the refilling port and insert a dropper containing fresh 3.33M potassium chloride (KCl) solution. The filling solution should reach the bottom of the refilling port.
- The filling solution level must be maintained just below the refilling port and higher than the pH buffer or sample level during calibration and measurement. This creates a positive head pressure forcing the filling solution to leak into pH buffer or sample through the junction and preventing the reverse.
- Bubbles may form and get trapped within the solution of the sensing tip or reference chamber during transportation. This can affect the operation of your pH electrode. To dislodge the bubbles, gently shake the electrode body.
- If the filling solution inside the reference chamber gets contaminated with sample or microbial growth or the reading is drifting, change the filling solution. Tilt the pH electrode, uncover the refilling port, and draw out the old solution using a dropper before refilling it with fresh 3.33M KCl solution.

Conditioning

Nowadays, combination and 3-in-1 pH electrodes are commonly available. Both types of pH electrodes consist of glass electrode and reference electrode built in one body, but the latter is integrated with temperature sensor for detecting the temperature of the solution being measured.

The glass electrode has a silver-based electrical wire suspended in a neutral solution with KCl contained inside a special glass. The surface of the glass bulb or membrane at the tip of the electrode must be hydrated to function properly. This can be accomplished by immersing the glass membrane in an aqueous solution, where a hydrated layer that is responsible for the pH response of the glass, is developed.

Another component of the pH electrode that must remain hydrated is the junction of the reference electrode. The junction is made of porous material such as ceramic or sintered polyethylene, which allows filling solution of the electrode to leak into the solution being measured. Keeping the reference junction hydrated will prevent precipitation of KCl from the filling solution which may clog it and cause erratic or slow electrode response.

- All pH electrodes come with white protective cap. A sponge wet with pure water is positioned at the bottom of the cap to keep the glass membrane and junction moist. If you find KCl salts formed on the junction or refilling port of your pH electrode, simply rinse off using clean water. This KCl creep from the filling solution is normal.
- A dry pH electrode will give inaccurate reading in pH measurement. Condition a dry pH electrode by soaking the glass membrane and junction in pH 7.00, 4.01 buffer, or tap water for at least 1 hour to regenerate the hydrated layer. Note: High salt solutions such as 3.33M KCl and the like are not recommended for conditioning our pH electrodes. After conditioning, rinse the pH electrode with clean water and proceed with calibration.

Never touch the glass membrane with fingers as oil or dirt may coat the glass and interfere with measurement.

Cleaning

A clean, hydrated glass membrane and free-flowing junction are necessary in performing an accurate measurement of pH. The choice of cleaning solution should effectively remove all contaminants based on sample tested without damaging your pH electrode.

1. If the pH electrode is liquid-filled, uncover the refilling port.
2. Clean the tip of your pH electrode using the appropriate cleaning solution. Make sure that the glass membrane and junction are both immersed in cleaning solution.
 - General samples—Soak the pH electrode in diluted detergent solution for 5 to 10 minutes, while moderately stirring the solution. A strong cleaning solution is needed for clogged junction, stains, and electrodes exhibiting slow response. Soak the pH electrode in cleaning solution 220 or 0.1M HCl for at least 1 hour.
 - Oily samples—Soak the pH electrode in warm, diluted detergent solution for 5 to 10 minutes, while moderately stirring the

solution. Alternatively, rinse the pH electrode with methanol or ethanol. Note: Alcohol is only applicable for glass-body electrodes. Never use organic solvents such as alcohol, acetone etc. to clean any plastic-body electrode as they may damage the body and shorten the life span. Use of organic solvents will void the electrode warranty.

- Protein-containing samples—soak the pH electrode in cleaning solution 250 for at least 1 hour.
3. Rinse the pH electrode with clean water.
 4. If the pH electrode is liquid-filled, draw out the old filling solution from the reference chamber and refill it with fresh 3.33M KCl (See Refilling).
 5. Condition the pH electrode (See Conditioning).

If calibration with fresh pH buffers failed repeatedly and cleaning failed to restore the performance, replace the pH electrode with a new one.

Storage

pH electrodes must be clean before they are stored for any length of time.

1. If the pH electrode is liquid-filled, cover the refilling port with the slider to prevent evaporation of filling solution.
2. Wash the protective cap with clean water to wet the sponge and remove KCl salts.
3. Insert the pH electrode into the protective cap with wet sponge. The water will not dissipate easily as the cap fit snugly on the electrode body. This environment is enough to keep the glass membrane and junction moist. It is not necessary to fill the cap with clean water and soak the pH electrode tip.

Short-term storage:

Between measurements, the pH electrode can be soaked in pH 7.00 buffer or clean water (e.g., tap, distilled or deionized).

Scan to download





ORP Powders



220

250

300



230

Cleaning Solutions





pH Buffer Solutions

| Code | Part No. | Description | Volume |
|---------|------------|---|--------|
| 100-2 | 3200043639 | pH 1.68 Buffer Solution at 25°C | 500ml |
| 100-4 | 3200043638 | pH 4.01 Buffer Solution at 25°C | 500ml |
| 101-SU | 3200738717 | pH Buffer Set (pH 4.01, 7.00, 10.01, 3.33M KCl) | - |
| 100-7U | 3200738711 | pH 7.00 Buffer Solution at 25°C | 500ml |
| 100-9 | 3200043636 | pH 9.18 Buffer Solution at 25°C | 500ml |
| 100-10U | 3200738712 | pH 10.01 Buffer Solution at 25°C | 500ml |
| 12.46 | - | Available on request | 500ml |

Conductivity Standard Solutions

| Code | Part No. | Description | Volume |
|--------|------------|---|--------|
| 100-21 | 3200738713 | 84 μ S/cm Conductivity Standard Solution | 500ml |
| 100-22 | 3200738714 | 1413 μ S/cm Conductivity Standard Solution | 500ml |
| 100-23 | 3200738715 | 12.88 mS/cm Conductivity Standard Solution | 500ml |
| 100-24 | 3200738716 | 111.8 mS/cm Conductivity Standard Solution | 500ml |
| 103-S | 3200738718 | Conductivity Standard Solution Set (84 μ S/cm, 1413 μ S/cm, 12.88 mS/cm, 111.8 mS/cm) | - |



ORP Powders

| Code | Part No. | Description | Volume |
|--------|------------|---|-----------------|
| 160-51 | 3200043618 | 89 mV at 25°C (for 250ml solution)   | 10 sachets/pack |
| 160-22 | 3200043617 | 258 mV at 25°C (for 250ml solution)   | 10 sachets/pack |

pH/ORP Electrode Filling Solutions

| Code | Part No. | Description | Volume |
|------|------------|-------------|--------|
| 300 | 3200043640 | 3.33M KCl | 250ml |

pH Electrode Cleaning Solutions

| Code | Part No. | Description | Volume |
|------|------------|---|--------------|
| 220 | 3014028653 | For removing inorganic residues from glass membrane and liquid junction   | 2 x 50ml |
| 230 | 3200530494 | For removing inorganic and organic residues from glass membrane (30ml Solution A & 100ml Solution B) | 30ml & 100ml |
| 250 | 3200366771 | For removing protein residues from glass membrane and liquid junction | 400ml |



Calcium Ion Electrode Solutions



Chloride Ion Electrode Solutions



Fluoride Ion Electrode Solutions




Potassium Ion Electrode Solutions




Ammonia Ion Electrode Solutions








Nitrate Ion Electrode Solutions



| Ion Standard Solutions | | | |
|------------------------|------------|---|--------|
| Code | Part No. | Description | Volume |
| 500-NH4-SH | 3200697171 | 1000 mg/L Ammonium Ion Standard Solution | 500ml |
| 500-NH4-SL | 3200697172 | 100 mg/L Ammonium Ion Standard Solution | 500ml |
| 500-CA-SH | 3200697175 | 1000 mg/L Calcium Ion Standard Solution | 500ml |
| 500-CA-SL | 3200697176 | 100 mg/L Calcium Ion Standard Solution | 500ml |
| 500-CL-SH | 3200697167 | 1000 mg/L Chloride Ion Standard Solution | 500ml |
| 500-CL-SL | 3200697168 | 100 mg/L Chloride Ion Standard Solution | 500ml |
| 500-F-SH | 3200697163 | 1000 mg/L Fluoride Ion Standard Solution | 500ml |
| 500-F-SL | 3200697164 | 100 mg/L Fluoride Ion Standard Solution | 500ml |
| 500-NO3-SH | 3200697179 | 1000 mg/L Nitrate Ion Standard Solution | 500ml |
| 500-NO3-SL | 3200697180 | 100 mg/L Nitrate Ion Standard Solution | 500ml |
| 500-K-SH | 3200697183 | 1000 mg/L Potassium Ion Standard Solution | 500ml |
| 500-K-SL | 3200697184 | 100 mg/L Potassium Ion Standard Solution | 500ml |

| Ionic Strength Adjustors | | | |
|--------------------------|------------|---|--------|
| Code | Part No. | Description | Volume |
| 500-NH3-ISA | 3200697174 | Ammonia Ionic Strength Adjustor  | 500ml |
| 500-CA-ISA | 3200697178 | Calcium Ionic Strength Adjustor | 500ml |
| 500-CL-ISA | 3200697170 | Chloride Ionic Strength Adjustor | 500ml |
| 500-F-TISAB | 3200697166 | Fluoride Ionic Strength Adjustor | 500ml |
| 500-NO3-ISA | 3200697182 | Nitrate Ionic Strength Adjustor | 500ml |
| 500-K-ISA | 3200697186 | Potassium Ionic Strength Adjustor | 500ml |

| Ion Selective Electrode Filling Solutions | | | |
|---|------------|--------------------------------------|--------|
| Code | Part No. | Description | Volume |
| 500-NH3-IFS | 3200697173 | Ammonia Electrode Filling Solution | 500ml |
| 500-CA-IFS | 3200697177 | Calcium Electrode Filling solution | 500ml |
| 500-CL-IFS | 3200697169 | Chloride Electrode Filling Solution | 500ml |
| 500-F-IFS | 3200697165 | Fluoride Electrode Filling Solution | 500ml |
| 500-NO3-IFS | 3200697181 | Nitrate Electrode Filling Solution | 500ml |
| 500-K-IFS | 3200697185 | Potassium Electrode Filling Solution | 500ml |

| Accessories | | |
|---|------------|--|
| Code | Part No. | Description |
| FA-70A | 3200644455 | Integrated Electrode Stand (Height: 338mm) for bench meter |
| FA-70S | 3200382557 | Adjustable, free-standing electrode stand (Height: 384 mm) |
| FA-70L | 3200382560 | Long, free-standing electrode stand (Height: 450-650mm) |
| DP-70S | 3200528474 | Electrode stand for 100 Series and D-70, ES-70, OM-70 Series handheld meters (Height: 400mm) |
|  | 3200373991 | Arm for electrode stand FA-70A, FA-70S, & FA-70L |
|  | 3200373961 | Electrode holders, 2pcs (for mounting electrode with round cap on electrode stand arm) |
|  | 3200382477 | Electrode protection caps, 3pcs (for 9615S-10D, 9618S-10D, 9681S-10D pH electrode) |
|  | 3200043508 | Electrode protection caps, 5pcs (for 9621-10D, 9625-10D, 9630-10D, 9631-10D, 9632-10D, 6367-10D, 6377-10D, 6252-10D, 6261-10C, 1066A-10C, 1076-10C, 2060-10T, 9300-10D, 9382-10D, 3552-10D pH electrode) |
|  | 3200382482 | Electrode protection cap for long electrode (for 9680S-10D, 9480-10C pH Electrode) |

With over 60 years of engineering excellence, HORIBA's diverse range of water quality analyzers and electrodes are ideal for everyday laboratory needs through to the most demanding of applications. Visit our website for a wealth of useful information and water quality measurement tips to help you obtain the best results in your work.



Benchtop Meters

Developed using extensive feedback from users, our new LAQUA meters deliver the best solution for water quality analysis. Our LAQUA website features an online 'Selection Guide' to enable you to find the perfect LAQUA meter and electrode for your need.



Handheld Meters

In the lab, in the field or anywhere you need it. LAQUA Handheld meters are designed for use with one hand and with an IP67 waterproof rating and shock-resistant casing. Meters can be used for long periods, even in dark places, making it ideal for field measurements in rivers and lakes.



Pocket Meters

Analyzing water quality is simplified when using our LAQUAtwin range of meters. Designed to produce accurate and reliable results. Anyone, anywhere, at any time can measure samples easily with a LAQUAtwin meter. See just how good they are at our website.



Application Notes

LAQUAtwin pocket meters offer quick and convenient alternative to analyze important parameters with high accuracy. Several application notes are available at (<http://goo.gl/znwE6j>) detailing the use of LAQUAtwin and the results achieved for the respective applications. Additional application notes will be added when available.

SUPPORT HORIBA CUSTOMER SUPPORT SYSTEM

HORIBA offers a variety of services to help you conform to quality standards and international guidelines such as GLP, GMP and ISO

Technical Support

Please contact us with any technical questions about our products.

www.horiba.com/wq/support

User Support

Our support website is available for registered customers and features:

- Data collection software
- Instruction manual downloads
- Measurement tips, etc.

www.horiba.co.jp/register

Validation Support

Please contact us with any questions or requirements for your validation procedure.

- Traceability certification*
- IQ/OQ/PQ support*
- SOP guidance
- FAQ

*Optional services



Please read the operation manual before using this product to assure safe and proper handling of the product.

- The contents of this catalog are subject to change without prior notice, and without any subsequent liability to this company.
- The color of the actual products may differ from the color pictured in this catalog due to printing limitations.
- It is strictly forbidden to copy the content of this catalog in part or in full.
- All brand names, product names and service names in this catalog are trademarks or registered trademarks of their respective companies.
- Windows is a registered trademark of Microsoft Corporation in the United States and other countries.

<http://www.horiba.com>

e-mail: labinfo@horiba.com

HORIBA Instruments Incorporation

9755 Research Drive
Irvine, CA 92618 U.S.A.
Tel: +1 (949) 250-4811
Fax: +1 (949) 250-0924



Distributed by: ABQ Industrial LP USA
Tel: +1 (281) 516-9292 / (888) 275-5772 eFax: +1 (866) 234-0451
Web: <https://www.abqindustrial.net> E-mail: info@abqindustrial.net

Brochure HEA-EN-03-2018B